



Client Lab Supply Requisition

PLEASE BE SURE TO COMPLETE THIS SECTION SO WE CAN CONTACT YOU WITH QUESTIONS REGARDING YOUR ORDER.

Client Account Name: _____

Address: _____

Requested By: _____

Today's Date: _____

Specific Delivery Location: _____

Phone # or E-mail: _____

Government regulations require NorDx to distribute lab supplies in reasonable proportion to the type and quantity of lab tests received.

| Quantity | |
|----------|--|
| | Needles / Sharps Containers |
| | 21G x 1 1/4" Eclipse (48 per box) |
| | 22G x 1 1/4" Eclipse (48 per box) |
| | Needle Holder-Single Use Eclipse (250 per bag) |
| | 5cc syringe (100 per box) |
| | 10cc syringe (100 per box) |
| | 2 Gal. sharps (ea.) |

| | |
|--|--|
| | Urine / Stool Collection |
| | 5 oz. urine container (ea.) |
| | 24 hr. container, plain (ea.) |
| | 24 hr. container, acetic acid (ea.) |
| | 24 hr. container, boric acid (ea.) |
| | 24 hr. container, HCl (ea.) |
| | IFOBT Personal Collection Kit - Mailer |
| | IFOBT - Vial Only |
| | Ova & parasite kits (ea.) |
| | 90 ml tox/pain urine container (ea.) |

| | |
|--|---|
| | Microbiology/Molecular Pathology |
| | Blood culture prep kit (ea.) |
| | Blood culture Adult Anaerobic (ea.) |
| | Blood culture Adult Aerobic (.ea.) |
| | Pediatric blood culture (ea.) |
| | Aerobic / Anaerobic culturette dual swab - blue cap (ea.) |
| | Minitip culture swabs - green cap (ea.) |
| | Plain Swab for Rapid Strep Screen (ea.) |
| | Universal transport (Flock/UTM) Herpes/Viral cult (ea.) |
| | Swab / Liquid Stuarts Dual BactiSwab-red cap (ea.) |
| | COVID / Flu / RSV / Pertussis with swab (20 per box) |
| | CT / GC / Vaginitis Aptima Multitest Swab (50 per box) |
| | CT / GC / Vaginitis Aptima Urine (ea.) |
| | CT / GC Aptima Unisex Swab (male urethral swab) (ea.) |
| | Stool culture kit (ea.) |

| | |
|--|---|
| | Forms and Labels |
| | Results form NorDx shadow logo |
| | 521 - Standard Assay Request (25 per pkg) |
| | 522 - Surg / Path / PAP Requisition form (25 per pkg) |
| | 523 - Micro Requisition form (25 per pkg) |
| | 525 - Allergy Requisition form (25 per pkg) |
| | 526 - HLA form (25 per pkg) |
| | Label: Patient Name / Date (500 per roll) |
| | Label: Zebra, medium (2500 per roll) |

| Other Items - Please specify below | |
|------------------------------------|--|
| | |
| | |
| | |

| Quantity | |
|----------|---|
| | Blood Collection Tubes |
| | 2 mL pediatric lavender (100 per pkg.) |
| | 2.7 mL pediatric light-blue 3.2% Citrate(100 per pkg) |
| | 3 mL green lithium heparin, PST (ea.) |
| | 4 mL gray potassium oxalate(ea.) |
| | 4 mL green lithium heparin (ea.) |
| | 4 mL lavender EDTA (100 per pkg.) |
| | 4 mL pediatric plain red (ea.) |
| | 5 mL gold SST (100 per pkg.) |
| | 6 mL plain red (ea.) |
| | 7 mL royal-blue EDTA (ea.) |
| | 7 mL royal-blue no additive (ea.) |
| | 8.5 mL yellow ACD (sol A) (ea.) |
| | 10 mL green sodium heparin (ea.) |
| | MaterniT 21 kits-streck black / tan tubes (ea.) |
| | 3mL tan (Pediatric Lead Testing) (ea.) |

| | |
|--|--|
| | Cytology and Histology Collection |
| | Thin prep vials (ea.) |
| | Thin prep broom (25 per pkg.) |
| | Thin prep spatula / brush (25 per pkg.) |
| | Cytobrush (100 per pkg.) |
| | Spray fixative (ea.) |
| | V.C.E. slides (72 per pkg.) |
| | Cardboard slide holder (ea.) |
| | 20 mL formalin vial (ea.) |
| | 40 mL formalin vial (ea.) |
| | 60 mL formalin vial (ea.) |
| | 90 mL formalin vial (ea.) |
| | 120 mL formalin vial (ea.) |

| | |
|--|---|
| | Miscellaneous Collection Supplies |
| | Alcohol wipes (Box 200 ea.) |
| | Gauze, 2x2 (200 Per sleeve) |
| | Hygienic towelettes (100 per pkg.) |
| | Tourniquets-Latex Free (ea.) |
| | Blood transfer device (ea.) |
| | C & S transfer device (ea.) |
| | Glucose tolerance drink - orange 100 gm (ea.) |
| | Glucose tolerance drink - orange 75 gm (ea.) |
| | Glucose tolerance drink - orange 50 gm (ea.) |
| | Transfer pipettes (500 per box) |
| | Transfer Tubes and Caps (100 per pkg.) |
| | Metal Free Aliquot tubes |
| | Amber Tubes |
| | Specimen bags (50 per pkg.) |
| | Specimen bags large 13X18 (100 per roll) |

NordxClientsupply@mainehealth.org
or fax order to 207-396-7848
Please allow 48-72 hours for delivery