

## **SUPPLY REQUISITION**

Office/Facility name \_ \_ Date \_\_\_  ${\sf Address}\,\_$ Phone\_ Ordered by \_

PHLEBOTOMY		QTY	QTY REQUISITIONS, FORMS AND PATIENT EDUCATION		
Tube gold 5 mL (SST)	(order as ea)		Laboratory supply requisition	(order as e	
Tube red 6 mL (plain)	(order as ea)		General lab requisition form (#443 or #435)	(100/pk	
Tube lavender 3 mL (EDTA)	(order as ea)		Anatomical pathology requisition form (#477)	(100/pk	
Tube light blue 2.7 mL (3.2% citrate)	(order as ea)		Advanced technology (ATL) requisition form (#480)	(100/pk	
Tube dark green 6 mL (sod heparin)	(order as ea)				
Tube light green 4.5 mL (lith heparin PST)	(order as ea)	├	Infectious disease (Microbiology) requisition form (#490)		
		<b>—</b>	Allergen test requisition form (#495)	(100/pl	
Tube royal blue 6 mL (trace element - Red stripe)	(order as ea)	_	Additional test requisition (add on)	(order as e	
Tube royal blue 6 mL MAYO (metal free EDTA - purple strip		_	HIV consent form (English)	(order as e	
Tube royal blue 6 mL MAYO (no additive aluminum)(limit 1		<u> </u>	HIV consent form (Spanish)	(order as e	
Tube tan 3 mL (EDTA)	(order as ea)		Cystic fibrosis brochure	(order as e	
Tube pink 6 mL (EDTA)	(order as ea)		Quad PT education booklet (English)	(order as e	
Tube gray 6 mL (sod fluoride)	(order as ea)	1	Quad PT education booklet (Spanish)	(order as e	
Tube yellow 6 mL ACD (sol A)	(order as ea)		Patient home collection instructions (English)	(order as e	
Tube yellow 6 mL ACD (sol B)	(order as ea)		Patient home collection instructions (Spanish)	(order as e	
Tube microtainer green (lith heparin)	(25/bag)		Orange Epic/Epic Care link lab order slips	(order as e	
Tube microtainer gold (SST)	(25/bag)	QTY	CYTOLOGY	(0.00.00	
Tube microtainer red (plain)	(25/bag)	Q I I	Pap test thin prep vial w/ cervical broom	(25/pk	
Tube microtainer lavender (EDTA)	(25/bag)	1—	Pap test thin prep vial w/ cervical brosh & spatula		
Tube QuantiFeron	(order as ea)	1—	,	(25/pk	
Tube urine myoglobin (mayo)	(order as ea)	├	Cytolyt for non GYN small mouth	(order as e	
Multi draw vacutainer needle 21 - 1¼"	(48/box)	ऻ—	Cytolyt for non GYN large mouth	(order as e	
			Green top slide fixative	(order as	
Multi draw vacutainer needle 22 - 1¼"	(48/box)		Fine needle aspirate (FNA) kit	(order as e	
Vacutainer holder	(each)		UroVysion kits	(order as	
Clear aliquot tube	(25/bag)	QTY	GLUCOLA		
Brown/Amber aliquot tube	(25/bag)	$ldsymbol{\Box}$	Glucola orange (50 gram)	(6/pl	
MAYO metal free aliquot tube	(25/bag)		Glucola orange (100 gram)	(6/pl	
MAYO clear aliquot tube	(100/bag)		Glucola fruit punch (100 gram)	(6/pl	
TOXICOLOGY			Glucola lemon lime (50 gram)	(6/pl	
Blue top drug screen cups	(order as ea)	QTY	SPECIMEN BAGS	(0/ pi	
Drug screen cup seals (numbered)	(order as ea)	ייע		(100 /-1	
Drug screen seals (unnumbered)	(roll)	1	White specimen bag (ambient)	(100/pl	
Polymedco Fit Kit	(order as ea)	1—	Green specimen bag (refrigerated)	(100/pl	
MICROBIOLOGY			Blue specimen bag (frozen)	(100/pl	
Standard tip Eswab	(order as ea)	1—	Red specimen bag (stat)	(100/pl	
Nasopharyngeal Eswab	(order as ea)	1	Purple specimen bag (priority)	(100/pl	
Blood culture bottle AFB	(order as ea)		Biohazard specimen bag (large)	(order as e	
Blood culture bottle aerobic	(order as ea)	QTY	MISCELLANEOUS		
Blood culture bottle anaerobic	(order as ea)	<u> </u>	Breath test kit	(5/bo	
Blood culture prep kit (Chloraprep)	(order as ea)		Specimen log book	(order as	
Blood culture device	(order as ea)	1	Pediatric lead kit	(order as	
			Pediatric lead kit (FILTER PAPER ONLY)	(order as	
Blood transfer device	(order as ea)		72 hour stool container	(2/pl	
Hologic aptima swab kit multitest/vaginal #301162 (pink)			Fetal fibronectin collection kit	(order as	
Hologic aptima swab kit unisex #301041 (purple)	(order as ea)	$\vdash$	Panorama prenatal screen	(order as	
Hologic aptima urine kit M/F #301040 (yellow)	(order as ea)	$\vdash$	·		
Pinworm paddle	(order as ea)	$\vdash$	Horizon collection kit	(order as e	
C & S, O & P stool kit	(order as ea)	$\vdash$	Horizon/Panorama combo kit	(order as	
STD multi-collection kit	(order as ea)		Cortisol saliva kit	(order as	
UTM viral transport tube	(order as ea)	QTY	HISTOLOGY		
URINE COLLECTION			Container with formalin 20 mL	(order as	
Sterile orange top urine container	(order as ea)		Container with formalin 60 mL	(order as	
Sterile drange top arme container	(order as ea)		Container with formalin 180 mL	(order as	
24 hour urine container - May Limit 10	(Jiuci as ca)		Container with formalin 480 mL	(order as	
24 hour urine container - Max Limit 10	(order as aa)		Container dry 180 mL	(order as	
Urine hat - max limit 25	(order as ea)		•		
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw)	(order as ea)	$\vdash$	Container dry 32 02	(order as a	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK)	(order as ea) (order as ea)		Container dry 83 oz		
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL)	(order as ea) (order as ea) (6/pkg)		Container dry 83 oz	(order as	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL)	(order as ea) (order as ea) (6/pkg) (6/pkg)		Container dry 83 oz Container dry 1 gallon/lids	(order as e	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea)		Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution	(order as e	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL)	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea) (order as ea)	QTY	Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution  LABELS	(order as e (order as e (order as e	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea)	QTY	Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution	(order as e (order as e (order as e	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea) (order as ea)	QTY	Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution  LABELS	(order as e (order as e (order as e (order as e (ro	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea) (order as ea) (order as ea)	QTY	Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution  LABELS  Specimen label - Epic care link 1" core 300778  Manuel specimen labels	(order as e (order as e (order as e (ro	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea) (order as ea) (order as ea) (order as ea)	QTY	Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution  LABELS  Specimen label - Epic care link 1" core 300778  Manuel specimen labels Patient Registry 3x1 - 1" core 300770	(order as e (order as e (order as e (rc (rc	
Urine hat - max limit 25 Urine C & S kit (Gray tube with transfer straw) Pediatric urine collection kit (PUCK) Acetic acid (15 mL) Acetic acid (25 mL) Boric acid 4 - 2.5 gram tablets Sodium carbonate 5 grams Sodium hydroxide 10 mL Urine yellow vacutainer 10 mL	(order as ea) (order as ea) (6/pkg) (6/pkg) (order as ea) (order as ea) (order as ea) (order as ea)	QTY	Container dry 83 oz Container dry 1 gallon/lids Flow Media/Hank solution  LABELS  Specimen label - Epic care link 1" core 300778  Manuel specimen labels	(order as e (order as e (order as e	

ALLOW 5 DAYS FOR DELIVERY OF SUPPLIES Spectrum Health Integrated Fleet Team (SHIFT) 616.486.5983